




喜诺
Genobio



Virusee[®]

SARS-CoV-2 Antigen Rapid Test (Colloidal Gold) (Saliva/Swab)

Product code : CoVSLFA-01 1 test/kit

CoVSLFA-20 20 tests/kit

Obtain result in 15 min



Overview

SARS-CoV-2, also known as Novel Coronavirus, causes acute pneumonia and is easily transmitted person to person and person to objects. After emerging, patients infected by SARS-CoV-2 have been successively discovered all over the world in a short time, and is still spreading rapidly. Nowadays, the diagnosis method widely used in clinical is real-time PCR. However, nucleic acid detection often gives false negative results due to sampling methods, virus titers, etc. An accurate and rapid detection is crucial in weathering the SARS-CoV-2.

Epidemiology

Source of infection: The novel coronavirus infected patients, asymptomatic patients can also be a source of infection

Transmission: Respiratory droplets and close contact transmission (Mainly), Aerosol transmission

Susceptible population: Generally susceptible

Product Characteristics


Name	SARS-CoV-2 Antigen Rapid Test (Colloidal Gold)
Method	Colloidal Gold
Specification	20 tests/kit, 1 test/kit
Sample type	Saliva, Nasopharyngeal Swab, Oropharyngeal Swab
Detection time	15 min
Stability	Stable for 12 months at 2-30 °C
Sensitivity	96.23%
Specificity	99.26%

Product Advantages


- **Flexible:** The sample type is optional between saliva and swab, convenient and meets different needing
- **Rapid:** Obtain result within 15 min
- **Simple:** Easy to use, without complicated operation
- **Intuitive result:** There is no need for calculation and instrument, visual reading result
- **Economic:** Product can be transported and stored at room temperature, reducing costs
- **Low risk:** Testing saliva or swab samples, reducing the risk of sampling process

Operation Procedure

1. Sample Collection



Saliva Sample



Nasopharyngeal Swab


Oropharyngeal Swab

Swab Sample


- **Saliva sample:** Spit at least 1 mL CLEAR saliva into the sample collection bag. The bag must not touch the mouth.
- **Swab sample:** Use sterilized swab to obtain a nasopharyngeal swab sample or oropharyngeal swab sample.

2. Sample Preparation


- **Saliva sample:**



- **Swab sample:**



3. Test Procedure



Add 3 drops to the specimen well (s).

Read the result at 15 min.



Innovation for Better Health

Tel: +86(22)60978609 ext. 8605/8606/8608/8609

Fax: +86(22)60978609-804

E-mail: marketing@era-bio.com / customer@era-bio.com

Web: www.era-bio.com



 **Genobio Pharmaceutical Co., Ltd.**

Add: 2-1,2,3-101, No.2 Zone, Binlv Industrial Park, Binhai
Tourism Area, Binhai New Area, Tianjin 300480, P. R. China

 **Tianjin Era Biology Technology Co., Ltd.**

Add: 1601, Building C2, TEDA MSD-C, No. 79, 1st Avenue, Tianjin
Economic-Technological Development Area, Tianjin 300457,
P. R. China

Genobio is wholly-owned subsidiary of Era Biology Group