

CE-IVD Marked

all Allplex™

M. genitalium & drug-resistance Assays

Simultaneous detection and identification of *Mycoplasma genitalium* and mutations responsible for resistance of azithromycin & moxifloxacin using multiplex real-time PCR

Allplex™ MG & AziR Assay

• Azithromycin resistance

Allplex™ MG & MoxiR Assay

• Moxifloxacin resistance

World's 1st technology

Individual Ct value of multiple analytes in a single channel



HIGH SENSITIVITY & SPECIFICITY

Multiplex real-time PCR with high sensitivity and specificity by utilization of DPO™ and TOCE™ technologies





M. genitalium & drug-resistance Assays


Continuing treatment of MG without diagnosis of drug resistance in MG will increase prevalence of resistance and decrease cure rate

Mycoplasma genitalium (MG) infection is a major cause of urethritis in men and is associated with cervicitis, pelvic inflammatory disease, preterm birth, and spontaneous abortion in women.

Men	Women
15~25% Symptomatic NGU*	10~30% Clinical cervicitis
10~30% NCNGU**	Among STD clinic attendees, 40~75% are asymptomatic
40% Chronic NGU*	Complications : PID***, tubal factor infertility


*) Non-gonococcal Urethritis, **) Non-chlamydial, ***) Pelvic inflammatory disease

Treatment guidelines²⁾




Increasing macrolide and fluoroquinolone resistance in MG infection¹⁾


The prevalence of **macrolide resistant mutations in MG infection**



The prevalence of **fluoroquinolone resistant mutations in MG infection**



Resistance-guided therapy improves cure rate of MG infection up to 92%³⁾



Cure rate of MG infection (%)

MG diagnostic workflow using Seegene's MG & DR Assays



Single test allows detection and identification of MG and two drug resistances in the same day.

One time, Right result	Drug	
	1 st - line	2 nd - line
<input checked="" type="checkbox"/> MG >	S	S
<input checked="" type="checkbox"/> MG > <input checked="" type="checkbox"/> AziR >	R	S
<input checked="" type="checkbox"/> MG > <input checked="" type="checkbox"/> MoxiR >	S	R
<input checked="" type="checkbox"/> MG > <input checked="" type="checkbox"/> AziR > <input checked="" type="checkbox"/> MoxiR >	R	R

S) Susceptible, R) Resistant

Allplex™ MG & AziR Assay

REAL TIME PCR CE-IVD Marked

Simultaneous detection and identification of *Mycoplasma genitalium* & 6 mutations in 23 S rRNA gene related to azithromycin resistance

Analytes

- *M. genitalium* - A2058T - A2059G - A2059C
- A2058G - A2058C - A2059T - Internal control



Allplex™ MG & MoxiR Assay

REAL TIME PCR CE-IVD Marked

Simultaneous detection and identification of *Mycoplasma genitalium* & 6 mutations in parC gene related to moxifloxacin resistance

Analytes

- *M. genitalium* - G248A - G259A - G259T
- G248T - A247C - G259C - Internal control



Key features of Allplex™ MG & DR Assays


- Simultaneous detection and identification of MG and mutations responsible for resistances of azithromycin (6 mutations) and moxifloxacin (6 mutations)
- Informative data with identification of each mutation related to drug resistances
- Short TAT (3.5 hours) from extraction to the final results
- Automated data interpretation by Seegene Viewer

Specimens

- Genital swab
- Urine
- Liquid based cytology (e.g., ThinPrep® and Surepath™)

Seegene's powerful automation platform for complete test process

Convenient total workflow using Seegene's automation platform




- Interlocked with LIS
- Multi-Ct values in a single channel
- User-friendly automation system
- Automated and optimized data analysis software for multiplex assays


Compatible instrumentation (CE-IVD Marked)

- Automated Extraction & PCR Setup
Seegene NIMBUS
Seegene STARlet
- Automated Extraction
NucliSENS® easyMAG®
SEEPREP32™
- Real-time PCR
CFX96™ Dx

Result of Allplex™ MG & DR Assays



Multiple Ct values in a single channel




The result represents detection of MG in Quasar670 with Ct values of 32.74 and A2059G in Cal Red 610 with Ct value of 34.34 (Graph 2). Allplex™ MG & AziR Assay and Allplex™ MG & MoxiR Assay are able to provide the information of MG and mutations related to azithromycin and moxifloxacin resistances by analyzing individual Ct value using Seegene Viewer.

Accurate discrimination of MG-DR mutations

Allplex™ MG & AziR and Allplex™ MG & MoxiR Assays are designed to detect and differentiate azithromycin and moxifloxacin resistance-associated mutations of MG as well as MG infection in the same time.


Distribution of drug-resistance mutations in MG positive samples

54% Azithromycin resistance mutation




Allplex™ MG & AziR Assay

14% Moxifloxacin resistance mutation



Allplex™ MG & MoxiR Assay



A total of 96 MG positive samples were tested.

- Azithromycin resistance mutations were detected in 52 (54.3%) samples; A2059G; 35 (36.5%), A2058G; 10 (10.5%), A2058T; 7 (7.3%)
- Moxifloxacin resistance mutations were detected in 14 (14.6%) samples; G248T; 4 (4.1%), G259A; 4 (4.1%), G259T; 3 (3.0%), G248A; 2 (2.0%), A247C; 1 (1.0%)
- Multidrug-resistance mutations related to AziR and MoxiR revealed in 11%.

Superior performance


Allplex™ MG & AziR Assay		
Detection	No. of detection	Agreement
MG	97 / 97	100%
Resistance markers total	53 / 97	100%
A2058T	7 / 97	100%
A2058C	0	-
A2058G	10 / 97	100%
A2059T	0	-
A2059C	0	-
A2059G	36 / 97	100%

Allplex™ MG & MoxiR Assay		
Detection	No. of detection	Agreement
MG	96 / 96*	100%
Resistance markers total	14 / 96*	100%
G259T	4 / 14	100%
G259C	2 / 14	100%
G259A	1 / 14	100%
G248T	3 / 14	100%
G248A	0	-
A247C	4 / 14	100%

* 1 sample could not be amplified by sequencing

A total of 97 MG positive samples were tested.

- The results of Allplex™ MG & AziR Assay and MG & MoxiR assay showed 100% agreement with sequencing results.



Allplex™ MG & DR Assays

- Simultaneous detection and identification of MG and mutations responsible for resistances of azithromycin (6 mutations) and moxifloxacin (6 mutations)
- Informative data with identification of each mutation related drug resistance
- Short TAT (3.5 hours) from extraction to the final results
- Automated data interpretation by Seegene Viewer

References

- 1) Emerging Infectious Disease (EID) Journal. CDC. 2018;24(2):1-228.
- 2) Jensen JS, et al. 2016 European guideline on Mycoplasma genitalium infections. J Eur Acad Dermatol Venereol. 2016;30(10):1650-1656.
- 3) Read TRH, et al. Outcomes of Resistance-guided Sequential Treatment of Mycoplasma genitalium Infections: A Prospective Evaluation. Clin Infect Dis. 2019;68(4):554-560.

Ordering Information

Category	Product	Package Volume	Cat. No.
Allplex™	MG & AziR Assay	50 rxns	SD10169Y
		100 rxns*	SD10170X
		25 rxns*	SD10232Z
Allplex™	MG & MoxiR Assay	50 rxns	SD10233Y
		100 rxns*	SD10234X

Instrument	Type	Cat. No.
CFX96™ Dx	Real-time PCR _ Optical Reaction Module	1845097-IVD
	Real-time PCR _ Thermal Cycler	1841000-IVD
Seegene NIMBUS	Automated Extraction & PCR Setup	65415-03
Seegene STARlet	Automated Extraction & PCR Setup	67930-03
SEEPREP32™	Automated Extraction	SG71100
STARMag 96 X 4 Universal Cartridge Kit*	Nucleic Acids Extraction Reagent	744300.4.UC384
STARMag 96 ProPrep**	Nucleic Acids Extraction Reagent	EX00009P EX00009T

* For use with Seegene NIMBUS and Seegene STARlet. ** For use with SEEPREP32™

Not Available for Sale in the United States



Taewon Bldg., 91, Ogeum-ro, Songpa-gu, Seoul, Republic of Korea, 05548 / Tel : +82-2-2240-4000 / E-mail : info@seegene.com

www.seegene.com

BRAZIL

Belo Horizonte, Brazil
Tel : +55-31-25153003
E-mail : contato@seegenebrazil.com.br

CANADA

Toronto, Canada
Tel : +1-800-964-5680
E-mail : canada@seegene.com

GERMANY

Düsseldorf, Germany
Tel : +49-211-9943-4260
E-mail : sgg@seegene.com

MEXICO

México city, México
Tel : +52-55-8848-9646
E-mail : mexico@seegene.com

MIDDLE EAST

Dubai, UAE
Tel : +971-4-558-7110
E-mail : sgme@seegene.com

ITALY

Genova GE, Italy
Tel : +39-010-667-1796
E-mail : info@arrowdiagnostics.it

USA

California, USA
Tel : +1-925-448-8172
E-mail : info@seegenetech.com