

eazyplex® BloodScreen

When every minute counts — eazyplex® BloodScreen!
Help for therapy decision — in 20 min!




An early diagnosis and a prompt and appropriate clinical therapy of sepsis are crucial to increase the chance of survival. Antimicrobial resistance is a key factor determining clinical unresponsiveness to treatment and rapid progression in sepsis and septic shock.

With the eazyplex® BloodScreen tests as a Gram-positive and Gram-negative panel, you can detect the most common bacterial pathogens of sepsis and associated antibiotic resistances that strongly influence the outcome of the disease. You get the result without DNA extraction directly from positive blood culture in 20 minutes.

- Lyophilized tests
- 5 min sample pretreatment
- No DNA extraction
- Variable device selection
- Real-Time-Detection

eazyplex® BloodScreen GP A CE 0123 IVD	eazyplex® BloodScreen GN CE IVD
VanA	<i>E. coli</i>
VanB	<i>K. pneumoniae</i>
<i>Enterococcus faecalis</i>	<i>K. oxytoca</i>
<i>Enterococcus faecium</i>	<i>P. aeruginosa</i>
<i>S. pneumoniae</i>	CTX-M-1 group
<i>Streptococcus spec. I</i>	CTX-M-9 group
<i>Streptococcus spec. II</i>	<i>P. mirabilis</i>

Find out more and contact us at support@eazyplex.com



The eazyplex® test system is a freeze dried ready to use amplification reaction and can be stored directly at the workplace. During the amplification the result is shown in real-time.



AmplexDiagnostics GmbH
Werkstraße 2
83555 Gars Bahnhof
Germany
www.eazyplex.com

Kontakt/Contact

info@eazyplex.com
Tel.: +49 (0) 8073/91693-50
Fax: +49 (0) 8073/91693-55

