

# eazyplex<sup>®</sup> CSF direct

The first step in your meningitis diagnostics!  
Within 30 min you have the first important information!



- Lyophilized tests
- 5 min sample pretreatment
- No DNA extraction
- Variable device selection
- Real-Time-Detection
- **CE IVD**

Direct CSF\* diagnostics ! With the eazyplex<sup>®</sup> CSF direct panel!

\*Direct testing from cerebrospinal fluid!


An early diagnosis and a rapid beginning of the treatment of bacterial meningitis / meningo-encephalitis is very important. Because the lethality increases strongly with a decelerated therapy start.

For a rapid diagnosis within 30 minutes directly from CSF without DNA extraction and therefore a fast beginning of the treatment, the eazyplex<sup>®</sup> CSF direct panel.

Don't waste any time when it counts!

	eazyplex <sup>®</sup> CSF direct	eazyplex <sup>®</sup> CSF direct V	eazyplex <sup>®</sup> CSF direct M
HSV-1	x	x	
HSV-2	x	x	
VZV	x	x	
<i>N. meningitidis</i>	x		x
<i>S. pneumoniae</i>	x		x
<i>S. agalactiae</i>	x		x
<i>L. monocytogenes</i>	x		x
<i>H. influenzae</i>			x
<i>E. coli</i>			x

Find out more and contact us under [support@eazyplex.com](mailto:support@eazyplex.com)



The **eazyplex<sup>®</sup>** test system is a freeze dried ready to use amplification reaction and can be stored direct at the workplace. During the amplification the result is shown in Real-Time.



AmplexDiagnostics GmbH  
Werkstraße 2  
83555 Gars Bahnhof  
Germany  
[www.eazyplex.com](http://www.eazyplex.com)

Kontakt/Contact

[info@eazyplex.com](mailto:info@eazyplex.com)  
Tel.: +49 (0) 8073/91693-50  
Fax: +49 (0) 8073/91693-55

