

eazyplex[®] C. difficile

Step-by-step diagnostics or all in one! It's your decision!
Within 25 min you have all the details!




- Lyophilized tests
- 5 min sample pretreatment
- No DNA extraction
- Variable device selection
- Real-Time-Detection
- **CE IVD**

Avoid serious complications such as pseudomembranous colitis and sepsis and clarify promptly the reason for sudden onset of diarrhea, nausea and fever. For the benefit of your patients and to support your hygiene management.

	eazyplex [®] C. difficile basic	eazyplex [®] C. difficile classic	eazyplex [®] C. difficile complete
GDH		x	x
Toxin B	x	x	x
Toxin A			x
Binary Toxin			x
Inhibition control	x	x	x

Find out more and contact us under support@eazyplex.com



The **eazyplex[®]** test system is a freeze dried ready to use amplification reaction and can be stored direct at the workplace. During the amplification the result is shown in Real-Time.



AmplexDiagnostics GmbH
Werkstraße 2
83555 Gars Bahnhof
Germany
www.eazyplex.com

Kontakt/Contact

info@eazyplex.com
Tel.: +49 (0) 8073/91693-50
Fax: +49 (0) 8073/91693-55

