

eazyplex[®] C. difficile

Step-by-step diagnostics or all in one! It's your decision!
Within 25 min you have all the details!




- Lyophilized tests
- 5 min sample pretreatment
- No DNA extraction
- Variable device selection
- Real-Time-Detection
- **CE IVD**

Avoid serious complications such as pseudomembranous colitis and sepsis and clarify promptly the reason for sudden onset of diarrhea, nausea and fever. For the benefit of your patients and to support your hygiene management. Fast and reliable as a basic cost-effective information service, or as a comprehensive and complete solution.

	eazyplex [®] C. difficile basic	eazyplex [®] C. difficile complete
GDH		x
Toxin B	x	x
Toxin A		x
Binary Toxin		x
Inhibition control	x	x

The toxin-producing bacterium *Clostridioides difficile* is one of the main causes of antibiotic-associated diarrhoea and the cause of virtually all cases of pseudomembranous colitis. Although around 2% of all healthy adults carry *C. difficile*, many patients acquire the bacterium through hospital-acquired infections.

Find out more and contact us under support@eazyplex.com



The **eazyplex[®]** test system is a freeze dried ready to use amplification reaction and can be stored directly at the workplace. During the amplification the result is shown in Real-Time.



AmplexDiagnostics GmbH
Werkstraße 2
83555 Gars Bahnhof
Germany
www.eazyplex.com

Kontakt/Contact

info@eazyplex.com
Tel.: +49 (0) 8073/91693-50
Fax: +49 (0) 8073/91693-55

