

eazyplex[®] EHEC

When things have to go fast — eazyplex[®] EHEC!
You know more — after 25 min!



copyright original picture M. Rohde/Helmholtz Zentrum für Infektionsforschung-HZI, Braunschweig

Just work with your favoured parameter combination and differentiate with only one test run direct from the stool sample or as culture confirmation between:

- EHEC
- EPEC
- EIEC
- EAEC


You can read about the application in routine practice in the publication by Penzel, A. et al.



- Lyophilized tests
- 5 min sample pretreatment
- No DNA extraction
- Variable device selection
- Real-Time-Detection
-  

	EHEC complete	EHEC basic
Verotoxin 1	x	x
Verotoxin 2	x	x
Intimin	x	
Haemolysin	x	
EIEC/Shigella	x	
EAggEC	x	
Verotoxin 2f	x	

Find out more and contact us at support@eazyplex.com



The eazyplex[®] test system is a freeze dried ready to use amplification reaction and can be stored directly at the workplace. During the amplification the result is shown in Real-Time.



AmplexDiagnostics GmbH
Werkstraße 2
83555 Gars Bahnhof
Germany
www.eazyplex.com

Kontakt/Contact

info@eazyplex.com
Tel.: +49 (0) 8073/91693-50
Fax: +49 (0) 8073/91693-55

