

eazyplex[®] TyphiTyper

When you need to create clarity — eazyplex[®] TyphiTyper!
You get precision — after 20 min!



- Lyophilized tests
- 5 min sample pretreatment
- No DNA extraction
- Variable device selection
- Real-Time-Detection
- **CE IVD**

A total of 2500 serotypes of Salmonella are known and they are the main cause of bacterial gastroenteritis. Worldwide there are approximately 93.8 million cases of salmonellosis each year causing 155,000 deaths.


Most of the salmonellosis are provoked by the Salmonella serovars S. Typhimurium and S. Enteritidis, while typhoid is induced by S. Typhi and S. Paratyphi A, B or C.

Create clarity and simplify your Salmonella diagnostics — and most of all, save time!

With the eazyplex[®] TyphiTyper you can detect the parameters mentioned below:

- Species Salmonella enterica and the following serovars
 - Salmonella enterica Serovar Typhi,
 - Salmonella enterica Serovar Paratyphi A, B and C in the itemization
 - Salmonella enterica Serovar Choleraesuis

Find out more and contact us at support@eazyplex.com



The **eazyplex[®]** test system is a freeze dried ready to use amplification reaction and can be stored direct at the workplace. During the amplification the result is shown in Real-Time.



AmplexDiagnostics GmbH
Werkstraße 2
83555 Gars Bahnhof
Germany
www.eazyplex.com

Kontakt/Contact

info@eazyplex.com
Tel.: +49 (0) 8073/91693-50
Fax: +49 (0) 8073/91693-55

