

# eazyplex® VRE

Take care of your patients!

Prevent the spread and nosocomial infections of VRE!

- Lyophilized tests
- 5 min sample pretreatment
- No DNA extraction
- Variable device selection
- Real-Time-Detection
- **CE IVD**




With the use of the eazyplex® VRE system you support your hygiene management and keep your budget balanced.

You get the result about VanA or VanB resistance in 20 min, whether from rectal swabs, blood culture or as culture confirmation, without DNA-extraction, including all necessary controls.

And as a pure culture confirmation test eazyplex® VRE basic, maximum cost efficient.  
So easy, fast and reliable as all our tests of the eazyplex® system.

Find out more and contact us under [support@eazyplex.com](mailto:support@eazyplex.com)



The **eazyplex®** test system is a freeze dried ready to use amplification reaction and can be stored direct at the workplace. During the amplification the result is shown in Real-Time.



AmplexDiagnostics GmbH  
Werkstraße 2  
83555 Gars Bahnhof  
Germany  
[www.eazyplex.com](http://www.eazyplex.com)

Kontakt/Contact

[info@eazyplex.com](mailto:info@eazyplex.com)  
Tel.: +49 (0) 8073/91693-50  
Fax: +49 (0) 8073/91693-55

