



ERA BIOLOGY




FungiXpert[®]

Inflammatory Biomarker series

Product Code: CRPCLIA-01, HBPCCLIA-01, PCTCLIA-01, IL6CLIA-01, SAACCLIA-01

Applicable Department

- Rheumatology
- Infectious department
- Cardiology
- Oncology
- Gastroenterology
- Nephrology
- Pulmonology
- Dermatology
- Neurology
- Endocrinology



Overview

Inflammation is a natural response of the body to injury, infection, or tissue damage. The use of inflammatory biomarkers in clinical practice has increased in recent years, as they have been found to provide valuable information about disease risk, progression, and treatment response. Commonly used inflammatory biomarkers include C Reactive Protein (CRP), Interleukin-6 (IL-6), Heparin Binding Protein (HBP), Procalcitonin (PCT) and Serum Amyloid A (SAA).

The detection of inflammatory biomarkers can provide valuable insights into the underlying causes of inflammation and guide appropriate management strategies.


Product Characteristics

Basic Parameters	Heparin Binding Protein (HBP) Detection Kit (CLIA)	Interleukin-6 (IL-6) Detection Kit (CLIA)	C Reactive Protein (CRP) Detection Kit (CLIA)	Procalcitonin (PCT) Detection Kit (CLIA)	Serum Amyloid A (SAA) Detection Kit (CLIA)
Limit of Detection	5.0 ng/mL	2.0 pg/mL	0.2mg/L	0.05 ng/mL	1.0 mg/L
Sample Type	Plasma	Serum, Plasma			
Specification	12 Tests/kit				
Applicable instrument	Full-Automatic Chemiluminescence Immunoassay System (FACIS)				
Detection time	40 min				
Method	Chemiluminescence Immunoassay				
Stability	The kit is stable for 1 year at 2°C-8°C				

Product Advantages



Saving time is saving life! Obtain result within 60 min.
(Including time for sample pre-treatment inside the system)



All in one! Integrated reagents and consumables together, more convenient.



Fully automatically perform experiment process & data analysis.

Tel: +86(22)60978609 ext. 8605/8606/8608/8609
 Fax: +86(22)60978609-804
 E-mail: marketing@era-bio.com / customer@era-bio.com
 Web: www.era-bio.com



 **Genobio Pharmaceutical Co., Ltd.**

Add: 2-1,2,3-101, No.2 Zone, Binlv Industrial Park, Binhai Tourism Area, Binhai New Area, Tianjin 300480, P. R. China



Tianjin Era Biology Technology Co., Ltd.

Add: 8-1, East District, Datang Headquarters Base, Binhai-Zhongguancun Science Park, TEDA, Tianjin, China

Genobio is wholly-owned subsidiary of Era Biology Group