




Anyplex™

MTB/NTMe

Real-time Detection

CE-IVD Marked



- ✓ Simultaneous detection of MTB & NTM
- ✓ Various sample types available
- ✓ High specificity and sensitivity
- ✓ Whole process control for assay validity
- ✓ Application of UDG System to prevent carry-over contamination
- ✓ Quick and easy data interpretation by Seegene viewer software

○ **Specimens (Anyplex™ MTB/NTMe Real-time Detection)**

- Sputum*
- Cultured cell
- Bronchial washing
- Fresh tissue

*NIMBUS IVD and STARlet IVD have been validated for sputum only

○ **Compatible Instrumentation** (CE-IVD Marked)

- **Auto Extraction & PCR setup**
NIMBUS IVD (Hamilton)
STARlet IVD (Hamilton)
- **Real-time PCR**
CFX96™ (Bio-Rad)

○ **Analytes**

- MTB - Mycobacteria - Internal Control (IC)

Workflow-Automated diagnosis




MTB/NTMe Real-time Detection

REAL TIME CE-IVD Marked

Microbiological culture is a standard method used for diagnosis of active tuberculosis. However, culturing method has difficulties in distinguishing MTB from NTM due to their similar colony morphologies.

Anyplex™ MTB/NTMe Real-time Detection is a real-time PCR based product, which can differentiate *Mycobacterium tuberculosis* (MTB) with non-tuberculosis mycobacteria (NTM). Moreover, this product provides an Internal Control as a control for whole process. Therefore, it allows rapid and efficient screening of MTB and NTM to facilitate diagnosis and treatment of TB patients.

Result



Ordering Information

Not Available in U.S.A.

Product	Cat. No.	Package Volume
Anyplex™ MTB/NTMe Real-time Detection	TB7202Y	50 rxns
	TB7202X*	100 rxns

*For use with Microlab NIMBUS IVD & STARlet IVD only.



Seegene Inc.
Taewon Bldg. 91, Ogeum-ro, Songpa-gu, Seoul, 05548
Tel : +82-2-2240-4000 / Fax : +82-2-2240-4040
E-mail : info@seegene.com

Seegene Germany
Düsseldorf, Germany
Tel : +49-211-9943-4260
E-mail : eu@seegene.com

Seegene TECHNOLOGIES Inc.
California, USA
Tel : +1-925-332-5664
E-mail : usa@seegene.com

Seegene MIDDLE EAST
Dubai, UAE
Tel : +971-4-558-7110
E-mail : sgme@seegene.com

Seegene CANADA Inc.
Toronto, Canada
Tel : +1-800-964-5680
E-mail : canada@seegene.com

www.seegene.com