


eazyplex[®]

Pneumocystis jirovecii

Your fast diagnostic aid for immunocompromised patients!



- Lyophilized tests
- 5 min sample pretreatment
- No DNA extraction
- Variable device selection
- Real-Time-detection
- **CE IVD**


Direct diagnosis from BAL, BS and TS in 25 min.

With the eazyplex[®] Pneumocystis jirovecii you have a reliable support for your diagnostics!

P. jirovecii is the causative agent for Pneumocystis pneumonia (PCP), one of the most frequent and severe opportunistic infections in immunocompromised individuals. It is still under debate if PCP is a result from reactivation of latent childhood infection or from de novo infection. Of special interest for hygiene management is that a nosocomial patient-to-patient transmission was shown to be possible. To be able to treat a patient suffering from PCP as fast as possible, the eazyplex[®] Pneumocystis jirovecii can be performed out of bronchoalveolar lavage fluid (BAL), broncheal secretion (BS) and tracheal secretion (TS) without DNA extraction. The complete result will be shown in realtime and delivered as result report after 25 min.

Take the chance and step forward to the future of diagnostics with eazyplex[®].

Find out more and contact us under support@eazyplex.com



The **eazyplex[®]** test system is a freeze dried ready to use amplification reaction and can be stored direct at the workplace. During the amplification the result is shown in Real-Time.



AmplexDiagnostics GmbH
Werkstraße 2
83555 Gars Bahnhof
Germany
www.eazyplex.com

Kontakt/Contact

info@eazyplex.com
Tel.: +49 (0) 8073/91693-50
Fax: +49 (0) 8073/91693-55

AMPLEX
DIAGNOSTICS GMBH