


eazyplex® SuperBug CRE

Carbapenemases! Clear confirmation in 15 min!



- Freeze dried reactions
- 2 min sample pretreatment
- No DNA extraction
- Portable device
- Real-Time-detection
- **CE IVD**

Direct screening from rectal swab and urine in 20 minutes as well as confirmation from positive blood culture and culture plate in 15 minutes with eazyplex® CRE!

Due to excellent results the eazyplex® SuperBug CRE establishes its place in labs throughout Europe:

Tormo N. et al., ECCMID 2018

Rödel J., et al., EJCMI 2017

Cherkaou A. et al., SGM 2016

Hinić V., et al., JMM 2015

D'Inzeo T., et al., AMCLI 2015

Findlay J., et al., JAC 2015


García-Fernández S., et al., JAC 2015

See more under: www.eazyplex.com/publications

The detection of the following resistances is included:

- KPC
- NDM
- VIM
- OXA-48
- CTX-M-1 group
- CTX-M-9 group
- OXA-181

Find out more and contact us under support@eazyplex.com



The **eazyplex®** test system is a freeze dried ready to use amplification reaction and can be stored direct at the workplace. During the amplification the result is shown in Real-Time.

The purchase price of this product includes a non-transferable license under European Patents EP1020534, EP1873260, EP2045337, EP2287338, EP1275715, EP2253717, EP1327679, EP2031058 or their foreign counterparts, owned by Eiken Chemical Co., Ltd.



AmplexDiagnostics GmbH
Werkstraße 2
83555 Gars Bahnhof
Germany
www.eazyplex.com

Kontakt/Contact

info@eazyplex.com
Tel.: +49 (0) 80973/91693-50
Fax: +49 (0) 80973/91693-55

