

eazyplex[®] SuperBug expert

IMP, IMI, GES and GIM? Uncommon carbapenemases so far!




- Freeze dried reactions
- 2 min sample pretreatment
- No DNA extraction
- Portable device
- Real-Time-detection
- **CE IVD**

Confirmation panel for uncommon carbapenemases
in Enterobacteriaceae, *Pseudomonas spec.*, *Acinetobacter spec.*


Resistance to antibiotics in gram-negative pathogens is nothing new. Nevertheless, there is always a need for still rare targets in order to close a complete diagnosis.

With the test kit eazyplex[®] SuperBug expert you can quickly and reliably detect also the yet uncommon carbapenemases IMP, IMI, GES and GIM as fast culture confirmation in 20 minutes.



Source: Epidemiologisches Bulletin No 26 29th of June 2017, No. 28 12th of July 2018 and No 31 1st of August 2019.

Find out more and contact us under support@eazyplex.com



The **eazyplex[®]** test system is a freeze dried ready to use amplification reaction and can be stored direct at the workplace. During the amplification the result is shown in Real-Time.

The purchase price of this product includes a non-transferable license under European Patents EP1020534, EP1873260, EP2045337, EP2287338, EP1275715, EP2253717, EP1327679, EP2031058 or their foreign counterparts, owned by Eiken Chemical Co., Ltd.



AmplexDiagnostics GmbH
Werkstraße 2
83555 Gars Bahnhof
Germany
www.eazyplex.com

Kontakt/Contact

info@eazyplex.com
Tel.: +49 (0) 80973/91693-50
Fax: +49 (0) 80973/91693-55

