


eazyplex[®]

PneumoBug expert

All pneumonia is not the same!



- Lyophilized tests
- 5 min sample pretreatment
- No DNA extraction
- Variable device selection
- Real-Time-detection
- **CE0123 IVD**


Direct diagnosis from BAL, bronchial- and tracheal secretions in 25 min.
eazyplex[®] PneumoBug expert - the fast and reliable help for your diagnostics!

Fast and clear diagnosis of atypical pneumonia is a helpful tool for prompt decision-making on macrolide therapy. eazyplex[®] PneumoBug expert can be performed from bronchoalveolar lavage, bronchial secretions and tracheal secretions without DNA extraction. The result is displayed in real time and delivered as a result report after 25 minutes.

This is of particular advantage if the detection of pathogens such as *Chlamydia pneumoniae*, *Mycoplasma pneumoniae* and *Legionella pneumophila subsp. pneumophila / fraseri* was previously subject to very time-consuming and lab-intensive diagnostics.

Take the chance and step into the future of diagnostics with eazyplex[®].

Find out more and contact us under support@eazyplex.com



The **eazyplex[®]** test system is a freeze dried ready to use amplification reaction and can be stored direct at the workplace. During the amplification the result is shown in Real-Time.



AmplexDiagnostics GmbH
Werkstraße 2
83555 Gars Bahnhof
Germany
www.eazyplex.com

Kontakt/Contact

info@eazyplex.com
Tel.: +49 (0) 8073/91693-50
Fax: +49 (0) 8073/91693-55

AMPLEX
DIAGNOSTICS GMBH