

eazyplex[®] STD

STDs travel overnight. Stop it now!
Direct detection in 25 minutes!



- Lyophilized tests
- 5 min sample pretreatment
- No DNA extraction
- Variable device selection
- Real-Time-detection
- **CE0123 IVD**


Direct detection of STD causing pathogens from swab **in just 25 min!**

According to the WHO, one million people worldwide are infected with Sexually Transmitted Diseases every day.¹ This implies that timely diagnosis and therapy of this STD can not only prevent sequelae, but can also reduce the risk of transmission to sexual partners. Hence a quick and uncomplicated diagnosis is a good preventive tool to reduce the prevalence and further spread of STDs.

¹ [https://www.who.int/news-room/fact-sheets/detail/sexually-transmitted-infections-\(stis\)](https://www.who.int/news-room/fact-sheets/detail/sexually-transmitted-infections-(stis))

	STD complete	STD CT/NG	STD CT
<i>Chlamydia trachomatis</i>	x	x	x
<i>Neisseria gonorrhoeae</i>	x	x	
<i>Ureaplasma urealyticum</i>	x		
<i>Mycoplasma hominis</i>	x		
<i>Mycoplasma genitalium</i>	x		
<i>Treponema pallidum</i>	x		

Find out more and contact us under support@eazyplex.com



The **eazyplex[®]** test system is a freeze dried ready to use amplification reaction and can be stored direct at the workplace. During the amplification the result is shown in Real-Time.



Supported by:



on the basis of a decision
by the German Bundestag

AmplexDiagnostics GmbH
Werkstraße 2
83555 Gars Bahnhof
Germany
www.eazyplex.com

