


eazyplex[®] hv-K. pneumoniae

Your advantage for comprehensive patient management!



- Lyophilized tests
- 5 min sample pretreatment
- No DNA extraction
- Variable device selection
- Real-Time-detection
- **CE0123 IVD**

The most important pathogenicity factors for hypervirulence and hypermucoviscosity.

- hypervirulent K. pneumoniae (hvKp) can lead to dangerous infections in healthy young people
- hvKp are significantly more pathogenic than "classic" K. pneumoniae
- typical manifestations are liver abscesses, pneumonia, meningitis and multilocal infections
- effective surveillance of hvKp helps to reduce outbreaks and risks for patients and public health


Without DNA extraction, you receive the result within 20 minutes, directly from the colony. Based on the result, you can immediately determine the virulence score (Kleborate) of the pathogen.

You can read how helpful our eazyplex[®] hv-K. pneumoniae kit is in routine use in the publication by Jürgen Rödel et al. (2023)



Rödel et al

Find out more and contact us under support@eazyplex.com



The **eazyplex[®]** test system is a freeze dried ready to use amplification reaction and can be stored direct at the workplace. During the amplification the result is shown in Real-Time.



AmplexDiagnostics GmbH
Werkstraße 2
83555 Gars Bahnhof
Germany
www.eazyplex.com

Kontakt/Contact

info@eazyplex.com
Tel.: +49 (0) 80973/91693-50
Fax: +49 (0) 80973/91693-55

